de

Medi-Test Combi 5 N®

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en

Medi-Test Combi 5 N®

Test strips for rapid determination of blood, protein, nitrite, ascorbic acid, glucose and pH-value in urine

Use
Screening test for early detection and supervision of diabetes and for detecting infections and diseases of the kidneys and urinary tract.
Only for use by qualified personnel.

Instructions for use
Dip the test strip for approximately 1 second into the fresh urine. Draw it across the rim of the container to remove excess urine. After 30 to 60 seconds (leukocyte test field after 60–120 seconds) compare the test strip with the color scale. The best time for comparison is after 30 seconds. Color changes that take place after more than 2 minutes are of no significance. When tested the urine should not be older than 2 hours.

Principle
Blood: The detection is based on the pseudoperoxidative activity of hemoglobin and myoglobin, which catalyze the oxidation of an indicator by an organic hydroperoxide producing a green color.
Protein: The test is based on the "protein error" principle of indicators. The test zone is buffered to a constant pH value and changes color from yellow to greenish blue in the presence of albumin. Other proteins are indicated with less sensitivity.
Nitrite: Microorganisms, which are able to reduce nitrate to nitrite, are indicated indirectly by this test. The principle of Griess reagent is the basis of this test. The test paper contains an amine and a coupling component. A red colored azo compound is formed by diazotisation and subsequent coupling.
Ascorbic acid: The detection is based on the decoloration of Tillmans reagent. In the presence of ascorbic acid a color change takes place from blue to red.
Glucose: The detection is based on the glucoseoxidase-peroxidase-chromogen reaction. Apart from glucose, no other compound in urine is known to give a positive reaction.
pH: The test paper contains indicators which clearly change color between pH 5 and pH 9 (from orange to green to turquoise).

Evaluation – Sources of Error
Blood: The minimum sensitivity of the test strip is 5 erythrocytes/µL urine corresponding to approx. 0.015 mg hemoglobin/dL urine. Intact erythrocytes are indicated by fleckly discolorations of the test field. The color fields correspond to the following values:
0 (negative), ca. 5–10, ca. 50, ca. 250 Ery/µL resp. hemoglobin concentration out of ca. 10, ca. 50, ca. 250 Ery/µL
Normal concentrations of ascorbic acid (< 40 mg/dL) do not influence the test results. Falsely positive reactions can be produced by a residue of peroxide containing cleansing agents.
Protein: The minimum sensitivity of the test strip is 10 mg protein/dL urine. The color fields correspond to the following ranges of albumin concentrations:
negative, 30, 100, 500 mg/dL or negative, 0,3, 1,0, 5,0 g/L
Falsely positive results are possible in alkaline urine samples (pH > 9), after infusions with polyvinylpyrrolidone (blood substitute), after intake of medicaments containing quinine and also by disinfectant residues in the urine sampling vessel. The protein coloration may be masked by the presence of medical dyes (e.g. methylene blue) or beetroot pigments.
Nitrite: The test detects concentrations from 0.05 mg nitrite/dL urine. Every pink color indicates a bacterial infection of the urinary tract. The color intensity depends only on the nitrite concentration, but does not provide information about the extent of the infection. A negative result does not preclude an infection of the urinary tract, if bacteria which cannot produce nitrite are present. Falsely negative results can be produced by high doses of ascorbic acid, by antibiotics therapy and by very low nitrate concentrations in urine as the result of low nitrate diet or strong dilution (diuresis). Falsely positive results can be caused by the presence of diagnostic or therapeutic dyes in the urine.
Ascorbic acid: The color fields correspond to the following values:
0 (negative), 10(+), 20(++) mg/dL or 0 (negative), 0,6(+), 1,1(++) mmol/L
Only for information!
Glucose: Pathological glucose concentrations are indicated by a color change from green to bluish green. Yellow or greenish test fields should be considered negative or normal. The color fields correspond to the following ranges of glucose concentrations:
neg. (yellow), neg. or normal (greenish), 50, 150, 500, ≥ 1000 mg/dL or neg. (yellow), neg. or normal (greenish), 2,8, 8,3, 27,8, ≥ 55,5 mmol/L
The influence of ascorbic acid (vitamin C) has been largely eliminated. An inhibitory effect is produced by gentisic acid. Falsely positive reactions can be produced by a residue of peroxide containing cleansing agents.
pH: The pH value of fresh urine of healthy people varies between pH 5 and pH 6. The color scale gives a clear distinction of pH value between pH 5 and pH 9.

Quality Control in professional use
The performance of the test strips should be confirmed by use of positive and negative control solutions. Positive and negative controls should be analyzed once a day, whenever a new bottle of strips is opened, whenever a new lot of strips is started, and every 30 days to check storage conditions. Each laboratory should establish its own goals for adequate standards of performance, and should question handling and testing procedures if these standards are not met.

Reactive ingredients
(Quantity resp. activity/cm² at time of impregnation)

Blood: tetramethylbenzidine cumene hydroperoxide	31 µg 315 µg	Nitrite: sulfanilic acid quinoline derivative	95 µg 37 µg	Glucose: glucose oxidase peroxidase tetramethylbenzidine	7 U 1 U 96 µg
Protein: tetrabromophenol blue	10 µg	Ascorbic acid: 2,6-dichlorophenolindophenol	7 µg	pH: methyl red bromothymol blue	3 µg 10 µg

Directions
In any case, in order to establish a final diagnosis and prescribe an appropriate therapy, the results obtained with test strips should be verified with other medical results.
The effect of medicaments or their metabolic products on the test is not known in all cases. In case of doubt it is recommended not to take the medicaments and then repeat the test.
Only use well washed and clean vessels for urine collection. The presence of usual urine preservatives will not affect the test results.
Remove only as many test strips as are required, and reseal the container immediately after use. Do not touch the test pads. Avoid exposing the strips to sunlight and moisture. Store the container below + 30 °C in a dry place. The test strips are stable, when stored properly up to the date of expiry indicated.
The caps contain a non-poisonous and harmless desiccant. In case this desiccant is swallowed accidentally, then drink plenty of water.
Explanation of symbols can be found in the package insert.
Disposal: Please dispose all used test strips in accordance with your local laws and regulations.
Package units: Tubes of 50 and 100 test strips
Date of change: 05/2015

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