




Medline Industries Ltd
3rd Floor, Quayside Wilderspool Business
Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom

Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk



Technical Sheet

UROSID[®] Silicone Foley Catheter

Transurethral Pediatric Silicone Foley Catheter, with stylet, sterile

Reference(s)

Reference	Size (CH)	Length (cm)	Balloon volume (ml)	Colour	Pack factor
641006	06	31	1.5	Light green	10
641008	08	31	3-5	Blue	10
641010	10	31	3-5	Black	10

Description

The UROSID[®] Silicone Foley Catheter is a double lumen, transurethral bladder catheter made of 100% silicone. It is designed with a Nelaton tip matching the needs of most of your patients. All UROSID Silicone Foley Catheters are equipped with a continuous barium sulfate X-ray contrast strip as well as a barium sulfate tip, to make sure the Foley catheter is well visualized during potential X-Ray examination.

CH6 - 10 are especially designed for pediatric use, offering urine drainage for your smallest patients. UROSID Foley Catheters are color-coded for easy selection of the correct CH-size. UROSID Silicone Foley catheters are equipped with a balloon to secure the catheter inside the patients' bladder.

The balloon has a capacity of 1.5ml up to 3-5ml, depending on the CH-size. For blocking of the balloon, a blocking medium is necessary. We recommend Glyco-SID (mixture of 10% glycerine and sterile water) for all UROSID Silicone Foley catheter. Only block the balloon with the maximum capacity mentioned on the blocking funnel using a Luer syringe.

Our pediatric Silicone Foley catheters come with a pre-inserted stylet for easy and safe positioning of the catheter.

Indwelling time: up to 30 days.

Not to be used in cases of:

Injuries of the urethra or stenosis, acute inflammation of the urinary tract (urethritis, prostatitis, epididymitis), in case of known hypersensitivity to the materials used

Main materials

Component	Material
Catheter	Silicone
Catheter tip	Silicone
Balloon	Silicone
Connector	(PVC, PP)
Stylet	(PA)

Pack factor

Refer to reference section



Retention period: End of life and support period

Packaging	Double sterile peel pouch
Legal Manufacturer	Well Lead Medical Co., Ltd. C-4 Jinhu Industrial Estate, Hualong 511434 Panyu, Guangzhou PEOPLE'S REPUBLIC OF CHINA In compliance with EN ISO 13485
Regulation	In compliance with Medical Devices Regulation 2017/745 and UK MDR MD Class: IIb sterile •Please refer to complete labeling before using the product.
Other Information	Sterile Ethylene oxide (EO) Shelf Life 5 years Storage -15°C to +49°C, in a dry place, away from sunlight