

Technical Data and security sheet Kidney dish

Réf. N° : 16 200 16/ 16 200 19/ 16 200 25 / 16 200 28
TCE1

Created :13/05/2013
Update :
Index : a
p 1/1

- Device class : Class I
- Period : Temporary (rule 1)
- No. of record ANSM : FR/CA01/0909763V
- GMDN Code : 12 143

Applied rule (annex IX)

- Non - Invasive
- Non surgical
- Non – diagnostic action
- Non therapeutic action
- Non Medical action

General description	Container made of stainless steel bean-shaped for receipt of medical waste.				
Enclosed documentation	Yes	Photo			
Batch reference number	Yes	On the back of the device			
Connected device	No				
Single Use Device	No				
Sterile Device	No				
Stérilization	Yes	Dry heat, autoclave and ethylene oxide.			
Autoclavable	Yes	134°C			
Cleaning	Yes	Soapy water or detergent which does not attack satinless steel			
Unsuitable substances	Products striker satinless steel				
Product components	No				
Materials	Stainless steel AISI 304				
Contains Latex	No				
Labeling	Yes	On the packing, designation, Name,Reference ,CE			
Packing	Yes	In an individual bag			
Particular storage	No				
Marking on the Medical Device	Yes	Reference, material, CE, Comed, batch number			
Références + extension		16 200 16	16 200 19	16 200 25	16 200 28
Désignation		Kidney dish 16 ml	Kidney dish 20 ml	Kidney dish 25 ml	Kidney dish 28 ml
Dimensions Length x Width x H (mm)		166 x 82 x 21	214 x 121 x 24	254 x 141 x 33	293 x 154 x 33
Capacity (ml)		200	500	800	1000
Net Weight (g)+- 1%		65.7	112.9	130.2	190
Colors		Polished satinless	Polished satinless	Polished satinless	Polished satinless
Temperature resistance (C°)		- 20 à +134°C	- 20 à +134°C	- 20 à +134°C	- 20 à +134°C
Use		Device used for temporary storage of instruments or similar products and for lipids			
Use care		Wash after each use risk of cross contamination			
Recycling		Yes			
Recovery		In Metal containers			
DEEE/REACH		No/No			
Phthalates/PVC		No/No			
Security data		R phrases : None Phrases combination = None			

Editor : Laura HEINRICH
Function : Quality Assistant
Paraph : LH

COMED sarl
8 Rue Louise Michel
67200 STRASBOURG France
www.comed.fr