




Biocompact®



Creados para facilitar su utilización e integrarse en los espacios hospitalarios aportando **seguridad y diseño**.

Para facilitar la extracción de agujas (Vacutainer y Venojet) hemos realizado un cómodo y seguro **acceso ambidiestro**.

Su exclusivo y practico **sistema de apertura y cierre** facilita su utilización, pudiendo realizarse sin ningún esfuerzo **con la presión de un simple dedo**.

Con el fin de facilitar su **colocación y manipulación** hemos dotado a esta gama de productos con **diferentes sistemas de fijación** que se adaptan a cada necesidad o situación, aportando seguridad en la manipulación de los residuos.



Crées pour faciliter leur utilisation et pour s'intégrer aux espaces hospitaliers en apportant **sécurité et design**.

Pour faciliter l'extraction des aiguilles (Vacutainer et Venojet) nous avons réalisé un **accès ambidextre** commode et sûr.

Son **système d'ouverture et de fermeture** exclusif et pratique, facilite son utilisation, celles-ci pouvant se réaliser sans aucun effort **avec la pression d'un simple doigt**.

Afin de faciliter son **emplacement et sa manipulation** nous avons doté cette gamme de produits avec **différents systèmes de fixation** qui s'adaptent à chaque besoin ou situation, en apportant de la sécurité dans la manipulation des déchets.



Created to facilitate its use and integrate in the hospital environment providing **safety and design**.

To facilitate the extraction of needles (Vacutainer and Venojet) we have implemented a comfortable and safe access for both the right-handed and left-handed.

Its exclusive and practical **opening and closing system** makes it easy to use; it can be done effortlessly, **with the pressure of a single finger**.

With the purpose of facilitating **location and handling**, this range of products has been fitted with **different systems of fastening** which adapt to each requirement or situation, providing waste-handling safety.

Ziocompact®


**ENVASES PARA RESIDUOS DE RIESGO
BOITE A DECHETS A RISQUES
SHARPS DISPOSAL CONTAINERS**



BS7320:1990



0,6L



1L



1,8L



3L



5L



10L



17L



22L



10L.



5L.



1.8L.



1L.



0.6L.



**MONTAJE
ASSEMBLY**

Presionar y asegurarse del perfecto encaje de la tapa con el cuerpo.
Press down and check the lid is fitted perfectly with the body.



1. TAPA ORIGINAL
ORIGINAL LID



2. TAPA APERTURA TOTAL
FULL OPENING LID




3. TAPA MINICOMPACT
MINICOMPACT LID

**APERTURA Y CIERRE
OPENING AND CLOSURE**

1. TAPA ORIGINAL
ORIGINAL LID




Presionar todos los puntos de cierre.
Push down at all the closing points.



2. TAPA APERTURA TOTAL
FULL OPENING LID



3. TAPA MINICOMPACT
MINICOMPACT LID

**SOPORTE PARED/BARRA
WALL/BAR SUPPORT****SOPORTE MESA
TABLE SUPPORT****BANDEJA
TRAY**

		DIMENSIONES: DIMENSIONS:	CAPACIDAD TOTAL: TOTAL CAPACITY:	CAPACIDAD ÚTIL: USEFUL CAPACITY:	CARGA MÁXIMA: MAX. LOAD:
MINICOMPACT	0.6L.	9.5 x 9.5 x 13 cm	0.6 L.	0.464 L.	0.4 Kg.
	1L.	9.5 x 9.5 x 19.5 cm	1 L.	0.78 L.	0.65 Kg.
BIOCOMPACT	1.8L.	16 x 16 x 17 cm	1.8 L.	1.44 L.	1 Kg.
	3L.	16 x 16 x 23 cm	3 L.	2.4 L.	2 Kg.
	5L.	22 x 22 x 18.5 cm	5 L.	4 L.	3 Kg.
	10L.	22 x 22 x 30.5 cm	10 L.	8 L.	5 Kg.
	17L.	38 x 28 x 29 cm	17 L.	13.6 L.	7.5 Kg.
	22L.	38 x 28 x 37 cm	22 L.	17.6 L.	9 Kg.
BIOCOMPACT APERTURA TOTAL	5L.	22 x 22 x 18.5 cm	5 L.	4 L.	3 Kg.
	10L.	22 x 22 x 30.5 cm	10 L.	8 L.	5 Kg.
	17L.	38 x 28 x 29 cm	17 L.	13.6 L.	7.5 Kg.
	22L.	38 x 28 x 37 cm	22 L.	17.6 L.	9 Kg.