



Article Closing-Cap, yellow, Luer, sterile  
Art. No. FK 8502822  
Intended use Closing stopper (Luer) without injection site  
Class IIa  
Rule 2  
Number of pages 3  
Date of revision 21.10.2011  
Replaces revision from  
Issued by B. Heyne  
Function: Manager Regulatory Affairs

Signed

*B. Heyne*  
.....

Revision list

Date	No. of revision	Modification
21.10.2011	0	First release



Article Closing-Cap, yellow, Luer, sterile

Art. No. FK 8502822

## Parts List

Position No.	Description	Material
1	Pin fastening, Luer, yellow	ABS (Acrylnitrile-Butadiene-Styrene), Colourant
Packaging		
	Compound foil PA/PE 5/90 µm	PA/PE (Polyamide / Polyethylene)
	Paper Fresenius	Paper
	Carton 185 x 100 x 72 mm for 100 pcs.	Carton
	Sticker, size: 97 x 67,7 mm	Paper
	Carton 400 x 220 x 380 mm for 2000 pcs.	Carton
	Additional information	Paper
	Sticker, size: 104 x 148 mm	Paper
	Adhesive tape paper	Tape
	Indicator adhesive tape	Crepe paper



Article

Closing-Cap, yellow, Luer, sterile

Art. No. FK

8502822

## Drawing

