



dicoNEX
single use syringe
3-part luer-lock, sterile
TD-41

Attachment No. 3
point 2 pp. 2.3

Revision IV
2018.12.31

Page 1/1

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

INDICATION:

- 3-piece syringe for intravenous, intramuscular and hypodermic injections
- Indicated for drug withdrawal/application and body fluids samples collection

FEATURES:

- Transparent barrel made of polypropylene allowing perfect visualization of the content
- Green contrasting plunger
- Luer-lock nozzle
- Sealing in the form of a double ring on a plunger stopper situated on the top of plunger
- Smooth plunger sliding thanks to rubber sealing
- Limitation flange securing against accidental removal of the plunger
- Name of the syringe and manufacturer's logo placed on the barrel
- Clear, legible and permanent black graduation facilitating dosage, adequately:
 - every 0,1 ml for the volume of 3 ml
 - every 0,2 ml for the volume of 5 ml and 10 ml
 - every 0,5 ml for the volume of 20 ml
 - every 1 ml for the volume of 30 ml and 50/60 ml
- Single-use
- Non-pyrogenic, non-toxic
- Latex-free
- DEHP-free
- EO sterilized
- Shelf life: 5 years
- Packaging: 1 pc./blister-pack; 25/50/100 pcs/box



SIZE AND REFERENCE

CATALOGUE NUMBER	SIZE	INDIVIDUAL PACKAGING	INTERMEDIATE PACKAGING	BULK PACKAGING
003ML-3CZ-LL-BL	3 ml	blister-pack	100 pcs	30 x 100 pcs
005ML-3CZ-LL-BL	5 ml	blister-pack	100 pcs	24 x 100 pcs
010ML-3CZ-LL-BL	10 ml	blister-pack	100 pcs	16 x 100 pcs
020ML-3CZ-LL-BL	20 ml	blister-pack	50 pcs	16 x 50 pcs
030ML-3CZ-LL-BL	30 ml	blister-pack	50 pcs	16 x 50 pcs
050ML-3CZ-LL-BL	50/60 ml	blister-pack	25 pcs	16 x 25 pcs