

HSW NORM-JECT®

HSW SOFT-JECT®

HSW Insulin Syringes

HSW Special Syringes

HSW FINE-JECT®

HSW reusable syringes



HSW Insulin Syringes

Manufacturing Programme disposable products: HSW-Disposable Syringes and Needles

According to the latest DIN, EN, ISO and FDA standards • Quality management system certified according to ISO 13485:2003 / AC 2007 • CE-marked

3-part disposable syringes HSW SOFT-JECT® with fixed needles

- ▶ needle tube permanently attached in syringe
- ▶ no dead space
- ▶ safety stop
- ▶ available in sizes 0,5 ml and 1 ml (0,3 ml coming soon!) both with Insulin U-100 and TBC graduation
- ▶ suitable for dosage of very small volumes
- ▶ needles available in two different lengths (12 mm and 8 mm)
- ▶ single-blister packed: 100 pieces per inner box

**with fixed
needles**

Size	Pieces per Europallet	Pieces per Box	Pieces per Case	Item Code
3-part HSW SOFT-JECT no dead-space				
1 ml 100IU fixed needle 30G x 1/2	48.000	100	2.400	8300010170
1 ml 100IU fixed needle 30G x 5/16	48.000	100	2.400	8300010462
0,5 ml 100IU fixed needle 30G x 1/2	48.000	100	2.400	8300013525
0,5 ml 100IU fixed needle 30G x 5/16	48.000	100	2.400	8300010463

Size	Pieces per Europallet	Pieces per Box	Pieces per Case	Item Code
2-part syringe HSW NORM-JECT no dead-space				
1 ml 100IU without needle	54.000	100	1.800	4010.A00V0
1 ml 40IU without needle	54.000	100	1.800	4010.100V0
1 ml 100IU by-packed needle 27G x 1/2	54.000	100	1.800	40127A00V0
1 ml 40IU by-packed needle 27G x 1/2	54.000	100	1.800	40127100V0
3-part syringe HSW SOFT-JECT				
1 ml 100IU without needle	54.000	100	1.800	5010.A00V0
1 ml 40IU without needle	54.000	100	1.800	5010.100V0
1 ml 100IU by-packed needle 27G x 1/2	54.000	100	1.800	50127A00V0
1 ml 40IU by-packed needle 27G x 1/2	54.000	100	1.800	50127100V0
1 ml 100IU mounted needle 27G x 1/2	32.000	100	1.600	8300009348