



dicoSPIKE / dicoSPIKE Chemo
Withdrawal cannula
for multi-dose containers

TD-44

att.no.3 pt.2 subpt.2.3

edition III

date 2021-01-12

change of day: 2021-12-22

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02,06 fax. +48 32 370 38 08

INDICATION:

- Device for needle-free dissolving and / or withdrawing drugs from multi-dose packages (bottles, vials, etc.) - in both directions
- dicoSPIKE Chemo model is designed for medicines releasing dangerous aerosols such as cytostatics

FEATURES:

- Easy to use, ergonomic shape
- Body made of ABS (dicoSPIKE Chemo – made of copolymer) with knurled surface for a secure grip
- Sharp receiving spike, 20 mm long, enabling tight and safe connection with containers of different capacity
- Groove in the receiving spike, 10 mm long, guarantees complete emptying of the container
- Hermetic self-locking cap securing the connector against contamination
- Available in version with air filter and in version with air and particle filter
- Air filter protects against microbiological contamination and against harmful aerosols
- Particle filter is an effective barrier to protect drugs and liquids from contamination
- Female luer lock connector compatible with products with either male or female luer lock
- Free of latex, PVC and phthalates
- Non-pyrogenic
- Single-use
- EO sterilized
- Packaging: 1 pc./ paper-foil



SIZE AND REFERENCE

CATALOGUE NUMBER	CAP COLOUR	AIR FILTER	PARTICLE FILTER	PURPOSE	INTERMEDIATE PACKAGING	BULK PACKAGING
PPLZ	green	0,1 µm	-	standard medicines	100 pcs	10 x 100 pcs
PPLN	blue	0,1 µm	5 µm	standard medicines	100 pcs	10 x 100 pcs
PPLC	red	0,1 µm	5 µm	cytostatics	100 pcs	10 x 100 pcs